



ORDER INFORMATION ALVISION™ (continued)

		KAPLAN		KAPLAN JS	
					
		KAPLAN 4.0 PP	KAPLAN 4.5 PP	KAPLAN JS 3.5 PP	KAPLAN JS 4.0 PP
5 Pack	Length	100 cm	100 cm	100 cm	100 cm
	Ø				
	5 F	501110019301	501110019302	501110019304	501110019303
5 Pack	Length	110 cm	110 cm	110 cm	110 cm
	Ø				
	5 F	501110019309	501110019310	501110019312	501110019311
5 Pack	6 F	501110019313	501110019314	501110019316	501110019315

		KAPLAN MOD I		
				
		KAPLAN MOD I 4.0 PP	KAPLAN MOD I 4.5 PP	KAPLAN MOD I 5.0 PP
5 Pack	Length	100 cm	100 cm	100 cm
	Ø			
	5 F	501110019786	501110019787	501110019890
5 Pack	Length	110 cm	110 cm	110 cm
	Ø			
	5 F	501110019893	501110019894	501110019892
5 Pack	6 F	501110019793	501110019794	501110019792

* Check with your sales representative the availability

PP = PRO PERFUSION - CURVES WITH SIDEHOLES

Kaplan MOD I and II : one frontal side hole
Kaplan PP: one frontal side hole
Kaplan JS PP: two lateral side holes



ORDER INFORMATION ALVISION™

		CORONARY BYPASS	
			
		CORONARY BYPASS LEFT (LCB)	CORONARY BYPASS RIGHT (RCB)
5 Pack	Length	100 cm	100 cm
	Ø		
	4 F	501110019621	501110019622
5 Pack	Length	100 cm	100 cm
	Ø		
	4 F	501110019721	501110019722
5 Pack	Length	100 cm	100 cm
	Ø		
	4 F	501110019821	501110019822
5 Pack	7 F*	501110019921	501110019922

* Check with your sales representative the availability

TECHNICAL DATA

Category:	Interventional cardiology diagnostic catheters
Diameters:	4F, 5F, 6F and 7F
Shaft :	Multi-segments (with variable stiffness)
Proximal shaft material:	Braided Nylon
Braiding:	Stainless steel flat wires
Distal shaft material:	Radiopaque PEBA
Tip material:	Soft PEBA with radiopaque Barium Sulfate
Guidewire compatibility:	0.038"
Maximum pressure:	1200psi

Flow Rate							
Size	Catheter			Pressure			
	Curve	ID	Length	300psi	600psi	900psi	1200psi
4F	Selective	0.042"	100 cm	11ml/s	15ml/s	18ml/s	22ml/s
	Pigtail	0.042"	110 cm	11ml/s	15ml/s	19ml/s	22ml/s
	Selective	0.042"	125 cm	10ml/s	14ml/s	17ml/s	20ml/s
	Pigtail	0.042"	125 cm	10ml/s	14ml/s	17ml/s	21ml/s
5F	Selective	0.042"	80 cm	12ml/s	16ml/s	20ml/s	24ml/s
	Selective	0.047"	100 cm	14ml/s	19ml/s	24ml/s	28ml/s
	Pigtail	0.047"	110 cm	14ml/s	19ml/s	25ml/s	28ml/s
	Selective	0.047"	110 cm	13ml/s	18ml/s	24ml/s	26ml/s
6F	Selective	0.047"	125 cm	13ml/s	17ml/s	22ml/s	25ml/s
	Pigtail	0.047"	125 cm	13ml/s	18ml/s	23ml/s	26ml/s
	Selective	0.047"	80 cm	15ml/s	21ml/s	26ml/s	30ml/s
	Selective	0.057"	100 cm	23ml/s	33ml/s	39ml/s	41ml/s
7F	Pigtail	0.054"	110 cm	21ml/s	29ml/s	35ml/s	37ml/s
	Selective	0.057"	110 cm	23ml/s	32ml/s	38ml/s	41ml/s
	Selective	0.057"	125 cm	21ml/s	30ml/s	35ml/s	37ml/s
	Pigtail	0.054"	125 cm	19ml/s	27ml/s	33ml/s	36ml/s
7F	Selective	0.057"	80 cm	26ml/s	36ml/s	41ml/s	45ml/s
	Selective	0.070"	100 cm	38ml/s	45ml/s	45ml/s	45ml/s
	Pigtail	0.067"	110 cm	35ml/s	45ml/s	45ml/s	45ml/s
	Selective	0.070"	80 cm	43ml/s	45ml/s	45ml/s	45ml/s

www.alvimedica.com

Please contact your local sales representative for product availability.
The images shown within this document are of a purely illustrative nature and are not intended to depict reality.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions and instruction for use.

©2023, Alvimedica. All Rights Reserved.




Alvision™ Interventional Cardiology Diagnostic Catheter

Alvision™


INTERVENTIONAL CARDIOLOGY DIAGNOSTIC CATHETER


Coronary angiographic catheters designed for effective diagnostic procedures

- Outstanding torque response and catheter stability due to Flat Wire Braid design
- Large inner lumens allow for superior visualization
- Wide choice of shapes to master every challenge



ORDER INFORMATION ALVISION™ (continued)


JUDKINS LEFT							
							
JUDKINS LEFT (JL 3.0) JUDKINS LEFT (JL 3.5) JUDKINS LEFT (JL 4.0) JUDKINS LEFT (JL 4.5) JUDKINS LEFT (JL 5.0) JUDKINS LEFT (JL 6.0) JUDKINS LEFT ST (JL 4.0 ST)							
Length Ø		100 cm	100 cm	100 cm	100 cm	100 cm	100 cm
5 Pack	4 F	501110019601	501110019602	501110019603	501110019605	501110019606	501110019607
	5 F	501110019701	501110019702	501110019704	501110019705	501110019706	501110019707
	6 F	501110019801	501110019802	501110019803	501110019805	501110019806	501110019807
	7 F *	501110019901	501110019902	501110019903	501110019905	501110019906	501110019907
Length Ø				125 cm		125 cm	125 cm
5 Pack	4 F			501110019502		501110019503	501110019501
	5 F			501110019507		501110019508	501110019509
	6 F			501110019513		501110019514	501110019515


JUDKINS RIGHT						
						
JUDKINS RIGHT (JR 3.0) JUDKINS RIGHT (JR 3.5) JUDKINS RIGHT (JR 4.0) JUDKINS RIGHT (JR 4.5) JUDKINS RIGHT (JR 5.0) JUDKINS RIGHT (JR 6.0)						
Length Ø		100 cm	100 cm	100 cm	100 cm	100 cm
5 Pack	4 F	501110019608	501110019609	501110019610	501110019611	501110019612
	5 F	501110019708	501110019709	501110019711	501110019710	501110019712
	6 F	501110019808	501110019809	501110019810	501110019811	501110019812
	7 F *	501110019908	501110019909	501110019910	501110019911	501110019912
Length Ø				125 cm		125 cm
5 Pack	4 F			501110019504		501110019505
	5 F			501110019510		501110019511
	6 F			501110019516		501110019517

* Check with your sales representative the availability

PP = PRO PERFUSION - CURVES WITH SIDEHOLES ST = SHORT TIP

ORDER INFORMATION ALVISION™ (continued)


AMPLATZ LEFT			
			
AMPLATZ LEFT (AL 1) AMPLATZ LEFT (AL 2) AMPLATZ LEFT (AL 3)			
Length Ø		100 cm	100 cm
5 Pack	4 F	501110019614	501110019615
	5 F	501110019714	501110019715
	6 F	501110019814	501110019815
	7 F *	501110019914	501110019915

PIGTAIL			
			
PIGTAIL STRAIGHT (PIG) PIGTAIL 145 (PIG 145) PIGTAIL 155 (PIG 155)			
Length Ø		110 cm	110 cm
5 Pack	4 F	501110019637	501110019638
	5 F	501110019737	501110019738
	6 F	501110019837	501110019838
	7 F *	501110019937	501110019938
	Length Ø	125 cm	125 cm
5 Pack	4 F	501110019525	501110019519
	5 F	501110019526	501110019521
	6 F	501110019527	501110019523

* Check with your sales representative the availability

PP = PRO PERFUSION - CURVES WITH SIDEHOLES ST = SHORT TIP

ORDER INFORMATION ALVISION™ (continued)

SONES						
						
SONES I (SON I) SONES II (SON 2) SONES III (SON 3) SONES I (SON I) SONES II (SON 2) SONES III (SON 3)						
Length Ø		100 cm	100 cm	100 cm	80 cm	80 cm
5 Pack	4 F	501110019628	501110019629	501110019630	501110019631	501110019632
	5 F	501110019728	501110019729	501110019730	501110019731	501110019732
	6 F	501110019828	501110019829	501110019830	501110019831	501110019832
	7 F *	501110019928	501110019929	501110019930	501110019931	501110019932
	PP	SONES I PP (SON I PP)	SONES II PP (SON 2 PP)	SONES III PP (SON 3 PP)		
	4 F	501110019634	501110019635	501110019636		
	5 F	501110019734	501110019735	501110019736		
5 Pack	6 F	501110019834	501110019835	501110019836		
	7 F *	501110019934	501110019935	501110019936		

Q4		HOCKEY STICK		MOD. EXTRA BACK UP		PROGRES-SIVE RIGHT	
							
Q4		HOCKEY STICK		MOD. EXTRA BACK UP		PROGRESSIVE RIGHT	
Length Ø	100 cm	Length Ø	100 cm	Length Ø	100 cm	Length Ø	100 cm
5 Pack	4 F	501110019647	501110019643	501110019648	501110019644		
	5 F	501110019747	501110019743	501110019748	501110019744		
	6 F	501110019847	501110019843	501110019848	501110019844		
	7 F*	501110019947	501110019943	501110019948	501110019944		

* Check with your sales representative the availability

PP = PRO PERFUSION - CURVES WITH SIDEHOLES ST = SHORT TIP